

Sirona Schick 33 Intraoral Sensor



Enhanced digital diagnostics

Sirona Dental's Schick 33 is enhancing digital dental diagnostics by combining the industry-leading theoretical resolution (33 line pairs/millimeter) and powerful enhancement software in one easy-to-use system.

For practices contemplating the switch to digital intraoral radiography, or those already experiencing the benefits, read more about what Schick by Sirona has to offer you. Our focus is clearly embodied in our promise of performance, protection and partnership. Through our dedication to these tenets, we maximize your control over patient care, help protect the long-term value of your investment, and help elevate the professional contributions of every member of your clinical team.

Performance

With an industry-leading theoretical resolution of 33 line pairs per millimeter, Schick 33 digital sensors capture and convey a full set of data. This yields a powerful, professional radiograph with the diagnostic accuracy your patients deserve.

The dynamic image enhancer lets you modify image sharpness, increasing your viewing capabilities by ensuring that you see more of the diagnostic data captured. You can view images from a variety of specific clinical perspectives: general dentistry, endodontics, periodontics, hygiene and restorative dentistry. Save these preferences and quickly share images with patients, colleagues and insurers in whatever format is preferred.

Patients come in multiple sizes, and so does Schick 33. Three sensor sizes (0, 1, and 2) give you the flexibility to get the right image for your treatment

situation, whether you're treating an adult or pediatric patient.

Easy-to-change, replaceable cables minimize downtime. And 3-foot, 6-foot, and 9-foot cable lengths can accommodate multiple operator environments.

Protection

To ensure durability, reliability and accuracy, Schick 33 digital sensors are assembled in New York by a team of dedicated craftspeople, many who have been with the company for more than 20 years.

These sensors are part of the Sirona comprehensive family of digital-dentistry solutions. With Schick 33, you ensure that your investment is protected when you decide to upgrade your patient experience through planned integration with CAD/CAM, 3D imaging and treatment centers.

Each of our products is backed by Sirona's and Patterson Dental's service-and-support program, including a large domestic call center, local tech support, in-office training, a comprehensive warranty, and a 100-percent-satisfaction guarantee.

Partnership

When you invest in Schick 33 digital radiography solutions, you become more than just a Schick by Sirona user. You take your place as a partner in our ongoing efforts to elevate the knowledge, performance, and status of every dental practice and every clinical professional. Your entire staff gains access to a full range of education, training and certification programs. You tap into Sirona's global network of technology leadership. ■