

ScanX[®] Duo

IMAGING | UTILITY ROOM | MERCHANDISE

INTRAORAL DIGITAL IMAGING

SEE MORE. DO MORE. **SPEND LESS.**

SEE MORE.

ScanX Duo captures high resolution, detail-rich images for improved diagnostics.

- Scan up to 22 line pairs per mm and see unsurpassed digital clarity.
- Enlarge, invert, or sharpen your radiographs to maximum diagnostic information and enrich patient discussions about treatment plans.
- Automatically store and file electronic images for future use or transmittal.

DO MORE.

ScanX Duo streamlines procedures and increases office productivity.

- Work more efficiently with no workflow disruptions.
- ScanX Duo can be used by 2 operatories simultaneously — a ScanX exclusive!
- Utilizing In-Line Erase, Phosphor Storage Plates (PSPs) are automatically scanned and erased in one continuous cycle- intraoral images and viewable in seconds.
- Improve patient comfort with thin, flexible, wireless, reusable phosphor plates in intraoral sizes 0, 1, 2, 3 and 4.

SPEND LESS.

ScanX Duo, both practical and cost-effective, will pay for itself in less than a year.

- Eliminate the darkroom, film, hazardous chemicals, and expensive ongoing processor maintenance costs.
- Works seamlessly with your existing X-ray equipment.
- Reuse inexpensive Phosphor Storage Plates thousands of times.



ScanX[®] Duo

SPECIFICATIONS:

Part Number:	D1000
Dimensions (HxWxD):	12"H x 10"W x 10.5"D (30.5cm x 25.4cm x 26.6cm)
Weight:	21 lbs. (9.5 Kg)
Electrical:	100 - 240V AC, 50/60 Hz
Resolution:	22 LP/mm (1100 dpi) Depending on scanning model/pixel size selected.
Theoretical Resolution:	40 LP/mm (2000 dpi)
Ports:	USB
Laser Classification:	Class I laser Product, Compliance with FDA 21 CFR 1040.10 and IEC 60825-1.
Imaging Medium:	Two processing slots accept Phosphor Storage Plates. Intraoral sizes 0, 1, 2, 3 and 4.
Accessories Included:	Power Cord, Computer Connector (USB), 2 each size 2 PSP Guides, 1 each size 0 and 1 PSP Guides, 8 each size 2 PSPs; 1 Box Barrier Envelope Size 2 (Qty. 300); 1 PSP Table Mat, 1 Sample ScanX Cleaning Sheet; 1 Sample PSP Cleaning Wipe; CD with Drivers, Utilities and Installation Manual.
Warranty:	ScanX Duo is warranted to be free from defects for a period of 24 months through an authorized Air Techniques dealer.

COMPUTER RECOMMENDATIONS:

CPU Speed:	2.8 GHz Pentium IV or Higher
Operating System:	Windows XP Professional with Service Pack 3 or later for an Intel 32-bit processor; Microsoft Windows XP Professional 64-bit Edition with Service Pack 2 or later for an Intel 64-bit extended (x64) processor; Microsoft Windows Vista Business, Enterprise, or Ultimate with Service Pack 1 or later for an Intel 32-bit or an Intel 64-bit extended (x64) processor; Microsoft Windows 7 Professional, Enterprise, or Ultimate or an Intel 32-bit or an Intel 64-bit extended (x64) processor.
RAM:	2 GB
ROM Drive:	48x speed
Hard Drive:	120 GB (2GB free for image acquisition, the rest for program files and image storage.)
Video Display Adapter:	256 MB RAM
USB Port:	2.0



Scan QR Code to learn more about ScanX Duo

www.airtechniques.com

Corporate Headquarters

1295 Walt Whitman Road | Melville, NY 11747-3062
Phone 800.247.8324 | Fax 888.247.8481

Western Facility

291 Bonnie Lane, Suite 101 | Corona, CA 92880-2804
Phone 800.247.8324 | Fax 951.898.7646

PN 9922-333 Rev. C

