

Bluephase® G4

LED curing light

Polyvision™
INSIDE 

With intelligent
curing assistant

Efficient
Esthetics


ivoclar
vivadent®
passion vision innovation

Intelligent LED light – carefree curing

Patented Polyvision technology is your personal assistant when it comes to curing dental materials.

This automatic assistance system helps you achieve reliable curing results. The curing light automatically detects movement of the handpiece during the light-curing procedure. It alerts the user of the improper operation by vibrating and automatically extends the exposure time by 10%, if necessary. If the movement may prevent the material from curing properly, the light will automatically interrupt the exposure cycle so that it can be repeated. Bluephase G4 comes with active anti-glare protection to prevent the light from being activated in open space.

Reliable curing performance for
all light-curing materials due to
polywave LED

Short curing times
starting from 10 seconds due to a high
intensity of 1,200 mW/cm²

Wide 10-mm light guide for
time-saving single-exposure
curing procedures

Additional Pre-Cure program
for the removal
of excess cement

Backlit display
for optimized hygiene

Homogeneous
light emission
performance

Light-curing assistant
Polyvision™

**Efficient
Esthetics**

Contactless battery charging

Part Number

691258BU Bluephase G4

ivoclarvivadent.com
EXPLORE SHOP LEARN

© 2018 Ivoclar Vivadent, Inc. Ivoclar Vivadent, Bluephase, Polywave and Polyvision are trademarks of Ivoclar Vivadent, Inc.

United States
175 Pineview Drive
Amherst, New York
14228
800-533-6825
Fax: 716-691-2285

ivoclar
vivadent
passion vision innovation