

T-Edge Autoclave

The only gravity autoclave
upgradable to pre/post vacuum
by software!

- Larger chamber capacity and shorter cycle times
- Convenient and easy loading with door opening to 110°
- Available in 110V and 220V
- 10" and 11" Chambers



Designed to meet all current
sterilization standards ANSI/AAMI ST55


Tuttnauer
Innovation · Legacy · Partnership



Features:

- Optimal chamber capacity
- Convenient and easy loading with door opening to 110°
- Assorted water filling and drainage methods (top, front or automatic)
- Range of automatic one-button programs
- User-friendly touchpad and control panel
- Advanced documentation and traceability options
- Remote Wi-Fi monitoring via smartphone or tablet

Printer includes sterilization data (including expiration date).



Technical Specifications

Autoclave Model	Chamber Size (In)	Volume	Overall Dimensions D x W x H (In)	Number of trays	Connections	Standard Cassette Capacity	Item Supplies	Power Supply	Power Consumption	Weight
T-Edge 10 220V	18" depth x 10" diameter	6 gal/ 23 liters	22.8"x18.9" x19.23"	5	1 ethernet port, 4 USB ports, WiFi, intergrated automatic water filling and drain connection	3 Full + 3 Half	Drain tube, tray holder	230V, Frq 50-60Hz	2000W	110 lbs.
T-Edge10 110V	18" depth x 10" diameter	6 gal/ 23 liters	22.8"x18.9" x19.23"	5	1 ethernet port, 4 USB ports, WiFi, intergrated automatic water filling and drain connection	3 Full + 3 Half	Drain tube, tray holder	120V, Frq 50-60Hz	1700W	110 lbs.
T-Edge11 110V	18" depth x 11" diameter	7 gal/ 27 liters	22.8"x18.9" x19.23"	5	1 ethernet port, 4 USB ports, WiFi, intergrated automatic water filling and drain connection	4	Drain tube, tray holder	120V, Frq 50-60Hz	1700W	115 lbs.

Sterilization Programs

Class S	Icon	Cycle Program	Load Type	Cycle Time with Drying (minutes)
		Unwrapped Instruments 270° F	Unwrapped Instruments	20 Minutes
		Wrapped Pouches 270° F	Wrapped Pouches	32 Minutes
		Unwrapped delicate 250° F	Unwrapped Delicate Instruments	45 Minutes
		Handpieces 270° F	Handpieces	45 Minutes

Class B	Icon	Cycle Program	Load Type	Cycle Time with Drying (minutes)
		Unwrapped Instruments 273° F	Unwrapped Instruments	20 Minutes
		Wrapped Pouches 273° F	Wrapped Pouches	39 Minutes
		Unwrapped Delicate 250° F	Unwrapped Delicate Instruments	35 Minutes
		Wrapped Delicate 250° F	Wrapped Delicate Instruments	60 Minutes

Sterisil® Ac+ System

Convenient, on demand autoclave water, made right in your practice

Designed for equipment needing distilled water

- Eliminates the need for purchasing and storing distilled water
- Produces up to 50 gallons of distilled quality water per day
- Output water is suitable for any instrument washer or autoclave.

Ultra Safe, Ultraviolet Equipped

- Class B UV Treatment - Reduces the elevated microbial load that inevitable occurs over time when storing water in tanks
- Helps ensure processed water is safe for unit waterlines

Built For Convenience

- Can be mounted under a sink or on the wall in a sterilization center - User serviceable without a dental technician
- Audible and visual maintenance notifications
- Water can be delivered via a faucet or the Sterisil* Wand (optional) depending on the practice needs.



STERISIL®



PN# STERISIL23



Tuttnauer

Innovation · Legacy · Partnership

Tuttnauer USA Co. Ltd.

25 Power Drive, Hauppauge, NY 11788

Tel: +800 624 5836, +631 737 4850, Fax: +631 737 0720

E-mail: info@tuttnauerUSA.com, www.tuttnauerUSA.com

Learn from our experts.

Visit our blog: <https://tuttnauerusa.com/blog/>

