



core3d  
centres®



# 3Shape TRIOS® 3 **NEW**

## Digital Impression Scanning

Core3dcentres® are an authorised distributor and trainer for 3Shape TRIOS® - the world's leading dental scanner and software supplier. The new TRIOS® 3 is the 3<sup>rd</sup> generation TRIOS® with its three in one digital impression solution. An intraoral scanner for fast and easy 3D colour impression taking; an integrated intraoral camera, so there's no need to purchase one separately, and shade measurement while you scan, for more accurate and predictable results. TRIOS® enables you to be more efficient, more precise and your patients more comfortable.

The 3Shape TRIOS® 3 comes in two options, the TRIOS® 3 Cart or the Pod Solution with the option for a handle or pen grip.

### Indications:

- Anterior and posterior single crowns
- Inlays
- Onlays
- Veneers
- Bridges
- Implant Abutments
- Temporaries
- Digital Models for Orthodontics
- Post and Core Restorations

### Advantages:

- Full arch scan in under 30 seconds
- Real colour scanning
- Extreme Mobility - The light and handy Pod solution is easy to share among treatment rooms or different locations
- Live 3D visualisation - View the digital impression being built on your PC screen during scanning
- Single hand operation - The newly designed scanner mount enables dentists to pick up and park the scanner with one hand
- Integrated training centre & remote support - Update your skills with latest videos and tutorials. Connect directly to TRIOS® support experts while sharing your screen
- Modern easy-to-use touch screen - Smart-touch screen provides full control without cumbersome trackball or mouse operation
- Unique motion sensor interface - Like a gaming-controller, use the hand-held to virtually rotate and turn the digital impression
- Flexible connectivity with Wi-Fi and Bluetooth - Move the cart anywhere in the clinic while staying online

# 3Shape TRIOS® 3

Digital Impression Scanning

## For the laboratory

### What to send:


To receive a restoration from Core3dcentres, the laboratory will be required to send either a digitally designed file from 3Shape Communicate™ with reference to the shade required and any special requirements or a upper and lower model with bite registration and clear instructions of what has to be manufactured.

## For the clinician

### What to send:

To receive a restoration from Core3dcentres, the Clinician will be required to send either a digital file from their 3Shape TRIOS® Digital Impression Scanner, using the 3Shape TRIOS inbox or 3Shape Communicate™ or a standard impression using traditional techniques with reference to the restoration required, the shade and any special requirements to their dental laboratory.

## Your next step...



Asia Pacific Benelux Canada Germany Japan Scandinavia South Europe Switzerland United Kingdom USA